

SOIL > SHEAR STRENGTH

# DIRECT SHEAR TESTING MACHINE (Large, Conventional)

Code : T220 (Large)



*will be added soon ...*

## TECHNICAL SPECIFICATIONS

*will be added soon ...*

## ORDERING INFORMATION

| Item Name                                | Item Code |
|------------------------------------------|-----------|
| DIRECT SHEAR TESTER (Conventional Model) | T220A10XU |
| LARGE DIRECT SHEAR TESTER - 300 mm       | T220A30XU |
| LARGE DIRECT SHEAR TESTER - 500 mm       | T220A50XU |

|                                               |           |
|-----------------------------------------------|-----------|
| SHEAR BOX ASSEMBLY - Ø 60 mm                  | T220P001H |
| SHEAR BOX ASSEMBLY - Ø 100 mm                 | T220P002H |
| SHEAR BOX ASSEMBLY - 60 x 60 mm               | T220P003H |
| SHEAR BOX ASSEMBLY - 100 x 100 mm             | T220P004H |
| LOADING TOP CAP - Ø 60 mm(for DIRECT SHEAR)   | T220P005H |
| POROUS DISCS - Ø 60 mm (for DIRECT SHEAR)     | T220P006H |
| POROUS STONES - Ø 60 mm (for DIRECT SHEAR)    | T220P007H |
| CUTTING EDGE - Ø 60 mm (for DIRECT SHEAR)     | T220P008H |
| EXTRUSION DOLLY - Ø 60 mm (for DIRECT SHEAR)  | T220P009H |
| LOADING TOP CAP - 60x60 mm(for DIRECT SHEAR)  | T220P010H |
| POROUS DISCS - 60x60 mm (for DIRECT SHEAR)    | T220P011H |
| POROUS STONES - 60x60 mm (for DIRECT SHEAR)   | T220P012H |
| CUTTING EDGE - 60x60 mm (for DIRECT SHEAR)    | T220P013H |
| EXTRUSION DOLLY - 60x60 mm (for DIRECT SHEAR) | T220P014H |