


# FALLING HEAD PERMEABILITY SET

Code : T250



- Used to determine the permeability of fine-grained soils

## **STANDARDS**

EN 17892-11

## **SUPPLIED WITH**

- Permeameter stand with 3 tubes
- Plated-Steel Falling-head permeability cell - Ø 100 mm (1 ea)
- Water Reservoir Tank

**ORDERING INFORMATION**

| Item Name                                   | Item Code |
|---------------------------------------------|-----------|
| MANOMETER STAND with 4 TUBES (Falling Head) | T250P001H |
| RESERVOIR TANK (Falling Head)               | T250P002H |
| PERMEABILITY CELL (Falling Head) - Ø 4"     | T250P003H |
| PERMEABILITY CELL (Falling Head) - Ø 6"     | T250P004H |
| MANOMETER TUBE - Ø 6 mm                     | T250P005H |
| MANOMETER TUBE - Ø 8 mm                     | T250P006H |
| MANOMETER TUBE - Ø 14 mm                    | T250P007H |
| MANOMETER TUBE - Ø 22 mm                    | T250P008H |
| FALLING HEAD PERMEABILITY SET               | T250X00XH |