

SOIL > MOISTURE DETERMINATION

# PORTABLE MOISTURE / TEMPERATURE METER (Microlance)

Code : T462



- Used to determine the moisture percentage and the temperature of the soil samples and the fine aggregates with maximum particle diameter of 10 mm. It measures upto 1000 mm depth.
- The range of the moisture content that can be detected by the equipment is between 0% (Dry condition) to 35% with sensitivity of (0.5%).

## **TECHNICAL SPECIFICATIONS**

- Depth: 1 m
- Moisture: 0-35% (0.5 %)
- Temperature: -20 / +60 (0.5°C)

**ORDERING INFORMATION**

| Item Name                                | Item Code |
|------------------------------------------|-----------|
| MICROLANCE MOISTURE / TEMPERATURE TESTER | T462X00XB |