

UNIVERSAL > PHYSICAL PROPERTIES

# INTERNAL PLASTIC PIPE PRESSURE TESTING MACHINE (Automatic)

Code : U012



*will be added soon ...*

## **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

## **ORDERING INFORMATION**

| Item Name                                     | Item Code |
|-----------------------------------------------|-----------|
| INTERNAL PLASTIC PIPE PRESSURE TESTER - 8 bar | U012A008U |
| PRESSURE CONTROLLING UNIT - 25 bar            | U012X025U |