

# VOLUMETRIC FLASK

Code : W003



- A volumetric flask in engineering testing laboratories is used for accurately preparing and diluting solutions to a specific volume. Its precision makes it ideal for measuring exact concentrations in chemical, environmental, or material testing, ensuring consistency and accuracy in experiments.

## ORDERING INFORMATION

| Item Name                  | Item Code |
|----------------------------|-----------|
| VOLUMETRIC FLASK - 5 ml    | W003G005H |
| VOLUMETRIC FLASK - 10 ml   | W003G010H |
| VOLUMETRIC FLASK - 1000 ml | W003G01KH |
| VOLUMETRIC FLASK - 25 ml   | W003G025H |

VOLUMETRIC FLASK - 2000 ml W003G02KH

VOLUMETRIC FLASK - 3000 ml W003G03KH

VOLUMETRIC FLASK - 4000 ml W003G04KH

VOLUMETRIC FLASK - 50 ml W003G050H

VOLUMETRIC FLASK - 5000 ml W003G05KH

VOLUMETRIC FLASK - 6000 ml W003G06KH

VOLUMETRIC FLASK - 100 ml W003G100H

VOLUMETRIC FLASK - 150 ml W003G150H

VOLUMETRIC FLASK - 250 ml W003G250H

VOLUMETRIC FLASK - 500 ml W003G500H

## OTHER PHOTOS

