

PYCNOMETER (Gay Lussac)

Code : W014



- A pycnometer (Gay Lussac type) in engineering testing laboratories is used to determine the specific gravity or density of liquids or powders. It consists of a precisely calibrated glass container, often with a stopper, ensuring accurate measurements by accounting for temperature variations during testing. This is essential for density and quality control in chemical, material, and environmental testing.

ORDERING INFORMATION

Item Name	Item Code
PYCNOMETER (Gay Lussac) - 5 ml	W014G005H
PYCNOMETER (Gay Lussac) - 10 ml	W014G010H
PYCNOMETER (Gay Lussac) - 1000 ml	W014G01KH

PYCNOMETER (Gay Lussac) - 25 ml	W014G025H
PYCNOMETER (Gay Lussac) - 2000 ml	W014G02KH
PYCNOMETER (Gay Lussac) - 50 ml	W014G050H
PYCNOMETER (Gay Lussac) - 100 ml	W014G100H
PYCNOMETER (Gay Lussac) - 250 ml	W014G250H
PYCNOMETER (Gay Lussac) - 500 ml	W014G500H

OTHER PHOTOS

