

PIPETTE

Code : W018



- A pipette in engineering testing laboratories is used to accurately transfer small volumes of liquid from one container to another. It is essential for precise measurements in chemical, biological, or environmental tests, ensuring consistency and accuracy in sample preparation and analysis.

ORDERING INFORMATION

Item Name	Item Code
PIPETTE - 1 ml	W018G001H
PIPETTE - 2 ml	W018G002H
PIPETTE - 5 ml	W018G005H
PIPETTE - 10 ml	W018G010H
PIPETTE - 25 ml	W018G025H