

CRUCIBLE TONG

Code : W022



- Crucible tongs in engineering testing laboratories are used to handle hot crucibles or containers, especially when they are heated to high temperatures. They provide a safe way to move and manipulate hot materials during processes like melting, burning, or other high-temperature chemical reactions.

ORDERING INFORMATION

Item Name	Item Code
CRUCIBLE TONG - 250 mm	W022M250H
CRUCIBLE TONG - 300 mm	W022M300H
CRUCIBLE TONG - 450 mm	W022M450H