

EVAPORATING DISH

Code : W024



- An evaporating dish in engineering testing laboratories is used for evaporating solvents from a solution, leaving behind the solid residue. It is typically made of porcelain or heat-resistant glass, designed to withstand high temperatures during the evaporation process, commonly used in chemical analyses or material testing.

ORDERING INFORMATION

Item Name	Item Code
EVAPORATING DISH - Ø 30 mm	W024C015H
EVAPORATING DISH - Ø 70 mm	W024C050H
EVAPORATING DISH - Ø 80 mm	W024C075H
EVAPORATING DISH - Ø 100 mm	W024C140H

EVAPORATING DISH - Ø 125 mm

W024C300H