

GLASS MORTAR and PESTLE

Code : W028



- A glass mortar and pestle in engineering testing laboratories are used for grinding and mixing materials, particularly for substances that may react with porcelain or other materials. The smooth surface of glass ensures no contamination and is ideal for preparing fine powders or combining chemicals in precise quantities for testing.

ORDERING INFORMATION

Item Name	Item Code
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H

GLASS MORTAR and PESTLE - Ø 180 mm

W028G180H

OTHER PHOTOS

