

HYDROMETER

Code : W154



- A hydrometer in engineering testing laboratories is used to measure the specific gravity or density of liquids. It is commonly used in soil testing, chemical analysis, and material characterization to determine the concentration of substances or to assess properties such as salinity or alcohol content.

ORDERING INFORMATION

Item Name	Item Code
HYDROMETER (0.700 - 0.800 g/ml)	W154G007H
HYDROMETER (0.800 - 0.900 g/ml)	W154G008H
HYDROMETER (0.900 - 1.000 g/ml)	W154G009H
HYDROMETER (1.000 - 1.100 g/ml)	W154G010H
HYDROMETER (1.100 - 1.200 g/ml)	W154G011H

HYDROMETER (1.200 - 1.300 g/ml)	W154G012H
HYDROMETER (1.300 - 1.400 g/ml)	W154G013H
HYDROMETER (1.400 - 1.500 g/ml)	W154G014H
HYDROMETER (1.500 - 1.600 g/ml)	W154G015H
HYDROMETER (1.600 - 1.700 g/ml)	W154G016H
HYDROMETER (1.700 - 1.800 g/ml)	W154G017H
HYDROMETER (1.800 - 1.900 g/ml)	W154G018H
HYDROMETER (1.900 - 2.000 g/ml)	W154G019H