

# VACUUM PUMP

Code : G702



Used to vacuum specimens, remove air voids and speed up deairing process.

## AVAILABLE MODELS

| Code         | Item                     |
|--------------|--------------------------|
| G 702 X051 U | VACUUM PUMP - 51 lt/min  |
| G 702 X100 U | VACUUM PUMP - 100 lt/min |
| G 702 X128 U | VACUUM PUMP - 128 lt/min |

## ORDERING INFORMATION

| Code         | Item                    |
|--------------|-------------------------|
| G 702 X051 U | VACUUM PUMP - 51 lt/min |

---

|              |                          |
|--------------|--------------------------|
| G 702 X100 U | VACUUM PUMP - 100 lt/min |
|--------------|--------------------------|

|              |                          |
|--------------|--------------------------|
| G 702 X128 U | VACUUM PUMP - 128 lt/min |
|--------------|--------------------------|

|              |                      |
|--------------|----------------------|
| G 702 P001 H | HOSE for VACUUM PUMP |
|--------------|----------------------|