

# LOADCELL

Code : G905



- Used for load measurements and calibration purposes
- Conforms to ASTM and EN Class 1
- Suitable for both laboratory and field use
- Supplied with durable carrying case

## **SUPPLIED WITH**

- Touch Screen Digital Indicator:
  - Unit Selection (kN / N / kgf / lbf)
  - Taring / Zeroing Option
  - Calibration (Password-Protected)
  - Data can be transferred to PC for reporting and printing in real-time
  - Battery and Mains operated

- Durable Carrying Case
- Calibration Certificate (Internationally traceable)

## **ORDERING INFORMATION**

Code	Name
G 905 X005 D	LOADCELL with DIGITAL INDICATOR - 5 kN
G 905 X01k D	LOADCELL with DIGITAL INDICATOR - 1000 kN
G 905 X025 D	LOADCELL with DIGITAL INDICATOR - 25 kN
G 905 X02K D	LOADCELL with DIGITAL INDICATOR - 2000 kN
G 905 X03K D	LOADCELL with DIGITAL INDICATOR - 3000 kN
G 905 X050 D	LOADCELL with DIGITAL INDICATOR - 50 kN
G 905 X250 D	LOADCELL with DIGITAL INDICATOR - 250 kN
G 905 X500 D	LOADCELL with DIGITAL INDICATOR - 500 kN

## **OTHER PHOTOS**

