

SOIL AUGER SET

Code : T305



- Used to take samples from the surface or near-surface soils. It is also used to extract disturbed samples from a soil bulk for boring logs
- Consists of a T handle, extension tube and an auger head
- The auger head comes in different diameters, to be chosen depending on the soil type and the desired amount of soil to be extracted.

STANDARDS

ASTM D420 • ASTM D1452 • AASHTO T86 • AASHTO T202

ORDERING INFORMATION

Item	Code
SOIL AUGER SET - Ø 2.5"	T 305 X063 H

SOIL AUGER SET - Ø 80 mm	T 305 X080 H
--------------------------	--------------

SOIL AUGER SET - Ø 100 mm	T 305 X100 H
---------------------------	--------------

SOIL AUGER SET - Ø 150 mm	T 305 X150 H
---------------------------	--------------

T-HANDLE for AUGER SET	T 305 P005 H
------------------------	--------------

EXTENSION ROD for AUGER SET - 1 m	T 305 P006 H
-----------------------------------	--------------

AUGER HEAD - Ø 2.5"	T 305 P001 H
---------------------	--------------

AUGER HEAD - Ø 80 mm	T 305 P002 H
----------------------	--------------

AUGER HEAD - Ø 100 mm	T 305 P003 H
-----------------------	--------------

AUGER HEAD - Ø 150 mm	T 305 P004 H
-----------------------	--------------